

Corona Test Prescription & Consultation

Release 2.4.4

RELEASE NOTE

Date	Author, department and contact details	For	Copy to	Comments
24/11/2020	Maxime Daive CovidPcrTestPrescription@smals.be	Click here to enter text.	Click here to enter text.	Click here to enter text.

TABLE OF CONTENTS

1.	RELEASE	3
2.	NEW FEATURES.....	3
2.1.	Addition of the list of laboratories for the “federal platform 2.0”	3
2.2.	No automatic form submission with scanner.....	3
2.3.	Check on sample ID	3

1. Release

RELEASE NUMBER	2.4.4	
TYPE		URGENT
APPLICATION	Corona Test Prescription & Consultation	
VERSION		
DATE ACCEPTATION	24/11/2020	
DATE PRODUCTION	25/11/2020	

2. New features

2.1. Addition of the list of laboratories for the “federal platform 2.0”

- Now, you can choose between “Federal platform 2.0” and “Other laboratory” when requesting the analysis of a sample
 - For “Federal platform” selection, display of the 8 supported laboratories and the “Federal Platform”
 - For “other laboratory” selection, display of the button to search through all labs
- The choice of the last selected laboratory is still saved to be displayed by default.

2.2. No automatic form submission with scanner

- Scanners automatically do a “carriage return” after the scan.
- In order to avoid uncontrolled submission, we have deactivated the sending of the form with the “Enter” key on the keyboard.

2.3. Check on sample ID

- If the entered sample format does not match the expected format, the warning message *“The format of the sample number does not seem to be correct”* is displayed.
- This message is just a warning to the user but does not block the submission of the form.
- This check is only done if a laboratory of “federal platform” is selected